

BRD73954

Chemical Properties


CAS No. : 1440209-96-0

Formula: C₁₆H₁₆N₂O₃

Molecular Weight: 284.31

Appearance: no data available

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 years



Biological Description

Description	BRD73954, an effective and specific HDAC inhibitor, which is with IC ₅₀ of 36 nM and 120 nM for HDAC6 and HDAC8, respectively.
Targets(IC ₅₀)	HDAC
In vitro	As compared to the other class I and II histone deacetylases (HDACs) (HDAC1, 12 μM; HDAC2, 9 μM; HDAC3, 23 μM; HDAC4, >33 μM; HDAC5, >33 μM; HDAC7, 13 μM; HDAC8, 120 nM; HDAC9, >33 μM), BRD73954 maintains excellent selectivity toward HDAC6 and HDAC8 with IC ₅₀ of 36 nM and 120 nM, respectively. Treatment with BRD73954 results in a robust increase in α-tubulin acetylation. [1]

Solubility Information

Solubility	Ethanol: 5.7 mg/mL (20 mM)),with gentle warming DMSO: 28.4 mg/mL (100 mM), (< 1 mg/ml refers to the product slightly soluble or insoluble)
------------	--

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	3.5173 mL	17.5864 mL	35.1729 mL
5 mM	0.7035 mL	3.5173 mL	7.0346 mL
10 mM	0.3517 mL	1.7586 mL	3.5173 mL
50 mM	0.0703 mL	0.3517 mL	0.7035 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Reference

Olson DE, et al. J Med Chem. 2013, 56(11):4816-4820.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only· Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:36 Washington Street,Wellesley Hills,MA 02481