

$^1\text{H}$  NMR (400 MHz,  $\text{DMSO-}d_6$ )  $\delta$  ppm 2.96 (dd,  $J=18.24, 9.76$  Hz, 1 H) 3.12 - 3.17 (m, 5 H) 3.24 - 3.31 (m, 1 H) 3.48 - 3.62 (m, 1 H) 3.62 - 3.71 (m, 1 H) 3.71 - 3.98 (m, 6 H) 4.40 - 4.66 (m, 1 H) 4.66 - 4.96 (m, 2 H) 6.65 (s, 1 H) 6.71 (s, 1 H) 6.74 (s, 1 H) 6.86 (s, 1 H) 9.20 (s, 1 H) 9.27 (s, 1 H)

